



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

Reagent for quantitative determination of Aspartate amino transferase activity (AST)
[EC 2.6.1.1] in human serum and plasma

AST GOT (IFCC)

REF LP80505	R1 4 x 30 mL	R2 1 x 30 mL
REF LP80605	R1 4 x 100 mL	R2 1 x 100 mL

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

support@biolabo.fr

Latest revision: www.biolabo.fr



Made In France

I: corresponds to significant modifications

INTENDED USE

This reagent is designated for professional use in laboratory (manual or automated method).

It allows the quantitative determination of aspartate amino transferase (AST) [EC 2.6.1.1] to screen its level in human serum and plasma.


GENERALITIES (1) (2)

AST is distributed in all body tissues, but greatest activity occurs in liver, heart, skeletal muscle and in erythrocytes. Minimal activity occurs in skin, kidney and pancreas. Although serum levels of both AST and ALT become elevated whenever diseases processes affecting liver cells integrity (viral hepatitis, liver necrosis and cirrhosis), an increased AST activity in serum or plasma appears in more than 97% of cases of myocardial infarction. AST levels (and occasionally ALT) are also elevated in progressive muscular dystrophy, pulmonary emboli, acute pancreatitis...

PRINCIPLE (4) (5)

Method developed by Karmen and Al and optimised by Henry and al. (according to modified IFCC recommendations).

Reaction scheme is as follows:



The decrease in absorbance proportional to AST activity in the specimen is measured at 340 nm.

Absence of P_iP allows a better stability of working reagent.

REAGENTS

R1	BUF ENZ AST	Buffer Enzymes
L-Aspartate		275 mmol/L
MDH		≥1000 UI/L
LDH		≥ 500 UI/L
EDTA		6 mmol/L
Tris Buffer		105 mmol/L
pH à 30°C		7,80 ± 0.1
Stabilizer		

R2	COENZ AST	Coenzyme
Tris Buffer		20 mmol/L
NADH		≤ 1,4 mmol/L
2-Oxoglutarate		80 mmol/L
Stabilizer		

According to 1272/2008 regulation (CLP), these reagents are not classified as dangerous

SAFETY CAUTIONS

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
 - Verify the integrity of the contents before use.
 - Waste disposal: Respect legislation in force in the country.
 - All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.
- Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENTS PREPARATION

Ready for use.

STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, reagents are stable when stored and used as described in the insert:

Unopened,

- Until the expiry date stated on the label of the Kit.

Once opened:

- Transfer requested quantity, well recap vials and store at 2-8°C,
- 2 separated reagents are stable at least 6 months without contamination
- Discard any cloudy reagent or if absorbance of mixed reagent (R1+R2) is < 1.000 at 340 nm.

SPECIMEN COLLECTION AND HANDLING (2)

Unhemolysed serum. Do not use heparinised plasma

AST is stable in serum or plasma for:

- 24 hours at room temperature
- 28 days at 2-8°C
- At least for 1 year at -20°C.

Adding pyridoxal phosphate (0.1 mM) improves stability at room temperature to 7 days.

LIMITS (3) (6)

LDH contained in reagent allows, during pre-incubation step, reduction of endogenous pyruvate which would positively interfere.

Likewise, oxaloacetate, product of the reaction, is carboxylated into pyruvate. This one will also be consumed by LDH contained in reagent and will not interfere with AST determination.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIAL REQUIRED BUT NOT PROVIDED

- Medical analysis laboratory equipment.
- Spectrophotometer or Biochemistry Clinical Analyzer

Manufacturer	Expiry date	In vitro diagnostic	Storage temperature	Dematerialized water	Biological risk
Product Reference	See Insert	Batch number	Store away from light	Sufficient for	Dilute with

QUALITY CONTROL

- **REF** 95010 EXATROL-N Level I
- **REF** 95011 EXATROL-P Level II
- External quality control program

It is recommended to control in the following cases:

- At least once a run
- At least once within 24 hours
- When changing vial of reagent
- After maintenance operations on the instrument

If control is out of range, apply following actions:

1. Prepare a fresh control serum and repeat the test
2. If control is still out of range, use a new vial of fresh calibrator
3. If control is still out of range, use a new vial of reagent and reassay

If control is still out of range, please contact BIOLABO technical support or your local Agent.

REFERENCE INTERVALS (1) (2)

	(IU/L) 37°C
New-born	39-117
Infant	23-94
Adult	13-31

Each laboratory should establish its own normal ranges for the population it serves.

PERFORMANCES

On KENZA 240TX, 37°C, 340 nm.

Linearity Range: between 9 and 500 IU/L

Detection limit: approx. 5 IU/L

Precision:

Within-run N = 20	Low level	Normal level	High level	Between run N = 20	Low level	Normal level	High level
Mean (IU/L)	23.6	45.9	162.8	Mean (IU/L)	23.8	45.7	170.5
S.D. IU/L	1.0	1.6	2.6	S.D. IU/L	1.3	2.1	4.1
C.V. %	4.2	3.5	1.6	C.V. %	5.4	4.6	2.4

Analytical sensitivity: approx. 0.0063 abs/min for 10 IU/L

Interferences:

Total bilirubin	Negative interference from 418 µmol/L
Direct bilirubin	No interference up to 420 µmol/L
Ascorbic acid	No interference up to 2500 mg/dL
Glucose	No interference up to 1060 mg/dL
Turbidity	Positive interference from 0.133 OD
Haemoglobin	Positive interference from 133 µmol/L

Other substances may interfere (see § Limits)

Comparison studies with commercially available reagent:

Realised on human specimens (n=100) between 10 and 330 IU/L

$$y = 0.9527x + 1.6243 \quad r = 0.9985$$

On-board stability: 2 separate reagents are stable 60 days

Calibration Frequency 14 days

Make a new calibration when changing reagent batch, if quality control results are found out of the established range and after maintenance operations.

CALIBRATION

- **REF** 95015 Multicalibrator traceable to *ERM-AD457/IFCC*

The calibration frequency depends on proper instrument functions and on the preservation of reagent.

PROCEDURE

Manual method:

Let stand reagents and specimens at room temperature

Pipette in 1cm pathlength thermostated cuvette	
Reagent 1	800 µL
Reagent 2	200 µL
Bring at 37°C, then add:	
Calibrator, Control or Specimen	100 µL
Mix. Start a timer. Record initial absorbance after 60 sec at 340 nm. Record the absorbance again every minutes during 180 sec.	
Measure absorbance change per minute (Δ Abs/min).	

- 1- Performances with manual procedure should be validated by user.
- 2- Kenza applications and other applications proposal are available on request.

CALCULATION

With Seric Multicalibrator:

$$\text{AST Activity} = \frac{(\Delta\text{Abs/min}) \text{ Specimen}}{(\Delta\text{Abs/min}) \text{ Calibrator}} \times \text{Calibrator Activity}$$

With Theoretical Factor:

$$\text{Activity (U/L)} = \Delta\text{Abs/min} \times \text{Factor}$$

$$\text{Factor} = \frac{\text{VR} \times 1000}{6.3 \times \text{VE} \times \text{P}}$$

With:

VR = Total reactional volume (mL)

VE = Specimen volume (mL)

6.3 = Molar extinction coefficient for NADH at 340nm

P = Path length (cm).

Example, with manual Procedure.

(Path length 1 cm, 37°C, 340 nm):

$$\text{IU/L} = (\Delta\text{Abs/min}) \times 1746$$

$$\mu\text{Kat/L} = \frac{\text{IU/L}}{60}$$

REFERENCES

- (1) *TIETZ N.W. Textbook of clinical chemistry, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 652-657*
- (2) *Clinical Guide to Laboratory Test, 4th Ed., N.W. TIETZ (2006) p. 154-159*
- (3) *YOUNG D.S., Effect of Drugs on Clinical laboratory Tests, 4th Ed. (1995) p. 3-68 to 3-79*
- (4) *HENRY R. J. and al., Am J clin Path (1960), 34, 381-398*
- (5) *IFCC Method for L-Aspartate aminotransferase. J Clin. Chem. Clin. Biochem. (1986), 24, p.497-510.*
- (6) *M. MATHIEU and col. SFBC. Comité de Standardisation. Recommandations pour la mesure de l'activité catalytique de l'Aspartate aminotransférase dans le sérum à 30°C. Ann. Biol. Clin. 1976. 34. 291-297*