

AQUATHERM

HEATING MATTRESS

GINEVRI
70
1954-2024

MADE IN ITALY



GINEVRI
neo-tech for life

A soft and warm nest for the little patient.

At GINEVRI, we understand the importance of making a newborn, even a premature one, feel as if they are in the comfort of the maternal womb. Upon birth, a baby requires various treatments designed to facilitate their growth phase and, most importantly, to prevent any trauma that could impact their physical state.

This is where our “Aquatherm” system comes into play. It provides an additional beneficial heat source for the little patient, making it an invaluable tool during this critical stage.

The Aquatherm System: Portable, Lightweight, and User-friendly

Our “Aquatherm” system is portable, lightweight, and easy to use. It includes:

- A Control Unit, equipped with feet for stability on any flat surface, as well as a support for placement on vertical surfaces.
- A Heating Mattress, which can be placed in a crib, on a care bench, a baby changing table, inside an incubator, or anywhere a heated flat surface is needed.

Optional: The Baby Bundle

As an optional addition, our “Baby Bundle” is a flexible and easily shapeable aid for the newborn to find their natural position. It’s a gentle and effective solution designed with the baby’s comfort in mind.

Versatility of the Aquatherm System

Thanks to its heat storage property, the “Aquatherm” system can be used in conjunction with the “Nido” baby crib for internal transport between different wards within the hospital.

TECHNICAL SPECIFICATION

Power Supply	230 Vac, 50-60Hz
Maximum Power	150 VA
Class	II B (93/42CEE)
Safety Class	I BF (EN60601-1)
Total weight	4.2 Kg

DIMENSIONS

Control Unit	19 x 15 x 19 cm
Heating Mattress	35 x 61 x 3 cm



CONTROL UNIT

Unit that proportionally servo-controls the rate of heat supplied to the mattress in order to reach and maintain the temperature set by the operator. It is provided with controls and alarms to ensure at the same time easiness of use and maximum reliability of the system.



HEATING MATTRESS

The mattress is the result of a new technology "Aqualgel", that offers all the advantages of a soft and impressed surface around the newborn. It is effective in minimizing the vibration and heat loss towards the surface and able to store the heat even when the control unit is not powered or in case of internal transport between wards.

BABY BUNDLE (optional)

a soft and flexible aid, which can be mounted around the baby to help him

Nido + Aquatherm configuration

ACCESSORIES

Baby Bundle



Nido Cradle



G/NEVRI

neo-tech for life

GINEVRI srl
Via Cancelliera, 25/b
00041 Albano Laziale (Rome) - Italy
Phone +39 06 9345 100
e-mail: mail@ginevri.com

Export Department
Phone +39 06 93459 331
e-mail: export1@ginevri.com

ginevri.com

Authorized Distributor

The specifications in this catalogue are indicative. The company GINEVRI reserves the right to make changes, without further notice, to the products described within this catalogue in order to improve reliability, function or design.

Quality System
ISO 9001:2015
ISO 13485:2016

