

KENZA ONE

Analyseur automatique de biochimie
Automated biochemistry analyser



Nouveau !
New!

Intelligent, compact et optimisé pour petits et moyens laboratoires
Smart, compact and optimized for small to medium laboratories.

Caractéristiques principales

- Toute méthode d'analyse
- Système Ouvert ou Fermé
- 100 tests par heure
- 23 positions réactifs (20 réfrigérées)
- Gestion STAT (Toutes les positions, à tout moment)
- 12 positions échantillons
- LED Sample Management
- Volume échantillon : 2-100 μ L
- Volume réactif : 3-350 μ L
- Volume réactionnel : 220-350 μ L
- 50 cuvettes longue-durée
- Démarrage rapide
- Bilan Digestif, Hépatique, Lipidique, Martial Rénal, Cardiaque, etc...
- Nouveau logiciel intuitif et interactif avec accès facilité aux informations

Main features

- All analysis methods
- Open or Closed System
- 100 tests per hour
- 23 reagent positions (20 refrigerated)
- STAT Management (At any position, at any time)
- 12 sample positions
- LED Sample Management
- Sample volume : 2-100 μ L
- Reagent volume : 3-350 μ L
- Reaction volume : 200- 250 μ L
- 50 long-life cuvettes
- Rapid start
- Analysis Panels : Digestive, Liver, Lipidic, Iron, Kidney, Cardiac, etc...
- New intuitive & interactive Software with easy and smart access to information



Gamme de réactifs en conteneurs code-barrés
Range of reagents in barcoded containers

Code	Produit / Product	
K1501	ACIDE URIQUE	URIC ACID
K1002	ALBUMINE	ALBUMIN
K1507	ALT / TGP	ALT / GPT
K1523	AMYLASE	AMYLASE
K1505	AST / TGO	AST / GOT
K1553	BILIRUBINE DIRECTE	DIRECT BILIRUBIN
K1443	BILIRUBINE TOTALE	TOTAL BILIRUBIN
K1004	CALCIUM	CALCIUM
K1005	CHLORURES	CHLORIDE
K1106	CHOLESTEROL	CHOLESTEROL
K1206	CHOLESTEROL-HDL	HDL-CHOLESTEROL
K1416	CHOLESTEROL-LDL	LDL-CHOLESTEROL
K1207	CK-NAC	CK-NAC
K1107	CREATININE	CREATININE
K1108	FER	IRON Direct
K1110	GAMMA GT	GAMMA GT
K1209	GLUCOSE	GLUCOSE
K1217	Isoenzyme CK-MB	CK-MB Isoenzyme
K1011	L.D.H.	L.D.H.
K1212	MAGNESIUM	MAGNESIUM
K1214	PHOSPHATASE ALCALINE	ALKALINE PHOSPHATASE
K1015	PHOSPHORE Inorganique	Inorganic PHOSPHORUS
K1016	PROTEINES TOTALES	TOTAL PROTEINS
K1519	TRIGLYCERIDES	TRIGLYCERIDES
K1532	UREE U.V.	UREA U.V.

Bientôt Disponible : Immunoturbidimétrie / Coming Soon: ImmunoAssays



BIOLABO
GROUP



BIOLABO S.A.S.
 Les Hautes Rives F-02160 Maizy (France)
 Tel. +33 (0) 323 251 550
 www.biolabo.fr sales@biolabo.fr support@biolabo.fr

Distributeur :
Distributor:

Photos non contractuelles / Non contractual photos

version: 08/11/2019